



Medical refrigerator / Medizinkühlschrank
MED288

Tel.: +49 40 600 094 680 E-Mail.: info@wesco-navy.de

Medical refrigerator MED288

Total capacity.....0,155 KW

=====

Marine Execution

Sealed unit compressor.

Volume : 288 Ltrs.


Consumption : 0,155 kW

Inside (WxDxH) : 530 x 500 x 1000 mm




Temperatures : +2 - +20 °C

- Door with lock.
- Door hinged left or right.
- Outside white painted steel plates.
- Inside: 5 medical drawers.
- Automatic defrosting, water evaporation, safety device against temperatures below zero, digital Temp.-display
- Alarm when temperatures are changing.
- Complete and ready for use with 180 cm cable and plug.

Side view



Refrigerator
from FRONT



* Power cord

This is a domestic refrigerator.
If you want to fasten it to the wall,
we recommend:

1. Use steel angle bars according to drawing.
2. The bars shall be fastened to the side walls
of the refrigerator and the ships wall.

Wir behalten uns vor, die technischen Spezifikationen im
Interesse der Weiterentwicklung zu verändern.
We reserve the right to modify the specifications
if necessary.

Medizinkühlschrank MED 288 Medical refrigerator		SCALE		FORMAT		DRAW.	AW	27.05.19
				A3				
		THIS DRAWING REMAINS THE PROPERTY OF WESCO Navy				REV.	-	-
Tel.: +49 40 600 094 680 Fax.: +49 40 536 75 01								
E-mail.: Info@wesco-navy.de		WESCO Navy Hamburg						